ANIOSPRAY SF IP sterile

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- *Formaldehyde free*
- *Biocidal product with broad spectrum and quick drying effect*
- *Activity in compliance with European standards*
- *LAL Testing ≤ 1 UI/ml*
- *Double-bagged packaged*
- *Gamma irradiated*

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CHARACTERISTICS

Rapid disinfection of surfaces and equipment in the Pharmaceutical Industry (Clean rooms).

COMPOSITION

Ethanol (CAS N°. 64-17-5 - 319 mg/g), glutaraldehyde (CAS N°. 111-30-8 - 0.53 mg/g), didecyl dimethyl ammonium chloride (CAS N°. 7173-51-5 - 0.86 mg/g), excipients.

MICROBIOLOGICAL PROPERTIES

Bactericidal in 15 min, 20°C (EN 1276, EN 13697).
Yeasticidal in 5 min, 20°C (EN 1650).
Fungicidal in 15 min, 20°C (EN 1650).
Active against Herpes virus and Rotavirus in 5 min, 20°C.
Activity on field surfaces and strains (bacteria, yeasts, moulds) according to EN 13697.
Other activities : see scientific file.

Delivered, for each batch, with Certificate of gamma irradiation, Certificate of analysis, LAL Testing analysis and Sterility report.

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PHYSICO-CHEMICAL DATA

- Density at 20°C: approx.
- pH of product at 20°C: 9 to 12
- Non-corrosive towards metals (non oxidising).

INSTRUCTIONS FOR USE

Used undiluted on previously cleaned surfaces.

Remove the first layer of packaging before entering into the clean room. Remove the second layer of packaging on first use.

Spray evenly over the surface to be treated respecting the indicated contact time for the required antimicrobial activity. Apply in sufficient amount (30 to 40 ml/m²) and ensure that the surfaces remain wet throughout.

For frequency of use, please refer to the hygiene plan in place.

PRECAUTIONS FOR USE

Dangerous. Follow the instructions for use (in accordance with Directive 99/45/EC and its adapted versions). Use the biocides with precaution. Before use, read the label and the product information.

Store between +5°C and +25°C.

Biocidal product designed for surface disinfection (Group 1 - PT2 - Directive 98/8/EC) - Reserved for professional use

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12 x 400ml spraycans gamma irradiated double-bagged: 1456.179