



Datasheet

STAXS® Reference	Description
IS-1005	DOTCH® Puru® Ball-point, black, double bagged in Al/LLDPE, 5 pc/inner bag, 1 inner bags/outer bag, 4 outer bags/box, sterile

The supplier is responsible for the information and data mentioned on the datasheet.
This version is the latest version received from the supplier.

DOTCH[®] Puru[®] Ball-point

STAXS[®] Reference : **IS-1005**

DOTCH[®] Puru[®] Ball-point, black, double bagged in Al/LLDPE, 5 pcs/inner bag, 1 inner bag/outer bag, 4 outer bags/box, sterile



Features & advantages

- Bic[®] ball-point
- Verifiable quality from a guaranteed source
- Continuous bioburden testing
- Packed in an aluminium bag and in an LLDPE cleanroom bag
- Double packed in ISO 5 cleanroom (ISO 14644-1)
- Validated sterility SAL 10⁻⁶
- Recommended for ISO 5 cleanrooms and higher

Physical characteristics

Type	Fountain pen	Colour	Black
		Line width	0,7mm

Technical data

Shelf life	Years	2*
------------	-------	----

* From date of manufacture

Revision date : 05/12/2017



SCIENTIFIC CREDIBILITY

DOTCH® Puru® Ball-point

STAXS® Reference : **IS-1005**

Packaging

Product	Ref. No.	Pieces/ inner bag	Inner bags/ outer bag	Outer bags/box	Box size	Box weight
DOTCH® Puru® Ball-point	IS-1005	5	1	4	33x24x4,5cm	

Notes:

- Packed in an aluminium inner bag and an LLDPE outer bag
- Box liner present in the box
- Minimum sales unit: 1 box
- Sterilization indicator on the inner label and on the box
- Production date, lot number and expiry date can be found on the inner label and on the box

Storage

- Store in dry, clean conditions in the original packaging and keep out of direct sunlight, high temperatures and solvent vapors
- Store between 5°C and 40°C

Certification

- Produced in conformity with ISO 9001:2008
- Bioburden and dosing (Periodic Dose Audit) checked periodically - ISO 11137-2
- Sterility validated according to ISO 11137-2:2012 AAMI - TR 27:2001 (VDmax²⁵)
- Certificate of irradiation (min. irradiation of 25kGy)
- Shelf life of 2 years from date of manufacture - ISO 11137-2 and ASTM F 1980-07

Revision date : 05/12/2017



Datasheet

STAXS® Reference	Description
IS-1006	DOTCH® Puru® Ball-point, blue, double bagged in Al/LLDPE, 5 pc/inner bag, 1 inner bags/outer bag, 4 outer bags/box, sterile

The supplier is responsible for the information and data mentioned on the datasheet.
This version is the latest version received from the supplier.

DOTCH[®] Puru[®] Ball-point

STAXS[®] Reference : **IS-1006**

DOTCH[®] Puru[®] Ball-point, blue, double bagged in Al/LLDPE, 5 pcs/inner bag, 1 inner bag/outer bag, 4 outer bags/box, sterile



Features & advantages

- Bic[®] ball-point
- Verifiable quality from a guaranteed source
- Continuous bioburden testing
- Packed in an aluminium bag and in an LLDPE cleanroom bag
- Double packed in ISO 5 cleanroom (ISO 14644-1)
- Validated sterility SAL 10⁻⁶
- Recommended for ISO 5 cleanrooms and higher

Physical characteristics

Type	Fountain pen	Colour	Blue
		Line width	0,7mm

Technical data

Shelf life	Years	2*
------------	-------	----

* From date of manufacture

Revision date : 05/12/2017



SCIENTIFIC CREDIBILITY

DOTCH® Puru® Ball-point

STAXS® Reference : **IS-1006**

Packaging

Product	Ref. No.	Pieces/ inner bag	Inner bags/ outer bag	Outer bags/box	Box size	Box weight
DOTCH® Puru® Ball-point	IS-1006	5	1	4	33x24x4,5cm	

Notes:

- Packed in an aluminium inner bag and an LLDPE outer bag
- Box liner present in the box
- Minimum sales unit: 1 box
- Sterilization indicator on the inner label and on the box
- Production date, lot number and expiry date can be found on the inner label and on the box

Storage

- Store in dry, clean conditions in the original packaging and keep out of direct sunlight, high temperatures and solvent vapors
- Store between 5°C and 40°C

Certification

- Produced in conformity with ISO 9001:2008
- Bioburden and dosing (Periodic Dose Audit) checked periodically - ISO 11137-2
- Sterility validated according to ISO 11137-2:2012 AAMI - TR 27:2001 (VDmax²⁵)
- Certificate of irradiation (min. irradiation of 25kGy)
- Shelf life of 2 years from date of manufacture - ISO 11137-2 and ASTM F 1980-07

Revision date : 05/12/2017