



# Datasheet

STAXS® Reference	Description
IS-1003	DOTCH® Puru® Marker, permanent ink, black, double bagged in Al/LLDPE, 5 pcs/inner bag, 1 inner bag/outer bag, 4 outer bags/box, sterile

The supplier is responsible for the information and data mentioned on the datasheet.  
This version is the latest version received from the supplier.

# DOTCH<sup>®</sup>

SCIENTIFIC CREDIBILITY

## DOTCH<sup>®</sup> Puru<sup>®</sup> Marker

STAXS<sup>®</sup> Reference : **IS-1003**

**DOTCH<sup>®</sup> Puru<sup>®</sup> Marker, permanent ink, black, double bagged in Al/LLDPE, 5 pcs/inner bag, 1 inner bag/outer bag, 4 outer bags/box, sterile**



### Features & advantages

- Permanent ink
- Verifiable quality from a guaranteed source
- Continuous bioburden testing
- Packed in an aluminium bag and in an LLDPE cleanroom bag
- Double packed in ISO 5 cleanroom (ISO 14644-1)
- Validated sterility SAL 10<sup>-6</sup>
- Recommended for ISO 5 cleanrooms and higher

### Physical characteristics

Ink composition	Ethanol Resin Colourants	Colour	Black
-----------------	--------------------------------	--------	-------

### Technical data

Shelf life	Years	1*
------------	-------	----

\* From date of manufacture

Revision date : 13/12/2017

PO Box 576 – 8440 AN - Heerenveen – The Netherlands  
T +31 (0)513 677 076 . F +31 (0)513 677 957 . Info-nl@staxs.eu

STAXS<sup>®</sup>

Page 1/2



SCIENTIFIC CREDIBILITY

## DOTCH® Puru® Marker

STAXS® Reference : **IS-1003**

### Packaging

Product	Ref. No.	Pieces/ inner bag	Inner bags/ outer bag	Outer bags/box	Box size	Box weight
DOTCH® Puru® Marker	IS-1003	5	1	4	33x24x4,5cm	

### Notes:

- Packed in an aluminium inner bag and an LLDPE outer bag
- Box liner present in the box
- Minimum sales unit: 1 box
- Sterilization indicator on the inner label and on the box
- Production date, lot number and expiry date can be found on the inner label and on the box

### Storage

- Store in dry, clean conditions in the original packaging and keep out of direct sunlight, high temperatures and solvent vapors
- Store between 5°C and 40°C

### Certification

- Produced in conformity with ISO 9001:2008
- Bioburden and dosing (Periodic Dose Audit) checked periodically - ISO 11137-2
- Sterility validated according to ISO 11137-2:2012 AAMI - TR 27:2001 (VDmax<sup>25</sup>)
- Certificate of irradiation (min. irradiation of 25kGy)
- Shelf life of 2 years from date of manufacture - ISO 11137-2 and ASTM F 1980-07

Revision date : 13/12/2017



# Datasheet

STAXS® Reference	Description
IS-1004	DOTCH® Puru® Marker, permanent ink, blue, double bagged in Al/LLDPE, 5 pcs/inner bag, 1 inner bag/outer bag, 4 outer bags/box, sterile

The supplier is responsible for the information and data mentioned on the datasheet.  
This version is the latest version received from the supplier.

# DOTCH<sup>®</sup>

SCIENTIFIC CREDIBILITY

## DOTCH<sup>®</sup> Puru<sup>®</sup> Marker

STAXS<sup>®</sup> Reference : **IS-1004**

**DOTCH<sup>®</sup> Puru<sup>®</sup> Marker, permanent ink, blue, double bagged in Al/LLDPE, 5 pcs/inner bag, 1 inner bag/outer bag, 4 outer bags/box, sterile**



### Features & advantages

- Permanent ink
- Verifiable quality from a guaranteed source
- Continuous bioburden testing
- Packed in an aluminium bag and in an LLDPE cleanroom bag
- Double packed in ISO 5 cleanroom (ISO 14644-1)
- Validated sterility SAL 10<sup>-6</sup>
- Recommended for ISO 5 cleanrooms and higher

### Physical characteristics

Ink composition	Ethanol Resin Colourants	Colour	Blue
-----------------	--------------------------------	--------	------

### Technical data

Shelf life	Years	1*
------------	-------	----

\* From date of manufacture

Revision date : 13/12/2017

PO Box 576 – 8440 AN - Heerenveen – The Netherlands  
T +31 (0)513 677 076 . F +31 (0)513 677 957 . Info-nl@staxs.eu

STAXS<sup>®</sup>

Page 1/2



SCIENTIFIC CREDIBILITY

## DOTCH® Puru® Marker

STAXS® Reference : **IS-1004**

### Packaging

Product	Ref. No.	Pieces/ inner bag	Inner bags/ outer bag	Outer bags/box	Box size	Box weight
DOTCH® Puru® Marker	IS-1004	5	1	4	33x24x4,5cm	

### Notes:

- Packed in an aluminium inner bag and an LLDPE outer bag
- Box liner present in the box
- Minimum sales unit: 1 box
- Sterilization indicator on the inner label and on the box
- Production date, lot number and expiry date can be found on the inner label and on the box

### Storage

- Store in dry, clean conditions in the original packaging and keep out of direct sunlight, high temperatures and solvent vapors
- Store between 5°C and 40°C

### Certification

- Produced in conformity with ISO 9001:2008
- Bioburden and dosing (Periodic Dose Audit) checked periodically - ISO 11137-2
- Sterility validated according to ISO 11137-2:2012 AAMI - TR 27:2001 (VDmax<sup>25</sup>)
- Certificate of irradiation (min. irradiation of 25kGy)
- Shelf life of 2 years from date of manufacture - ISO 11137-2 and ASTM F 1980-07

Revision date : 13/12/2017